

NIC Kolloquium

# An Introduction to the XNAT Data Management

December 10<sup>th</sup>, 2020

**Philipp T. Neukam**, Christian Löschen, Dirk K. Müller &  
the SFB940/INF Project

# Resources at the NIC homepage

<https://nic-tud.de/research/ressourcen/IT-services/>

## IT AND DATA SERVICES

### NIC Data-Management Manuals

- ↓ [NIC Data Management Overview](#)
- ↓ [Uploading DICOM files \(zip-file\) to the XNAT \(German\)](#)
- ↓ [Uploading non-DICOM files \(e.g. logfiles\) to the XNAT](#)

### Specific NIC Software Tools

- ↳ [Python script to parse the NIC-specific DICOM file paths -> nic\\_path.py](#)

### Naming Conventions at the NIC

- ↓ [Logfile Naming \(ideally\)](#)
- ↓ [MRI Sequenz \(Series\) Naming](#)

### Resources within the TUD VPN (restricted access)

- ↳ [Access to the VPN \(Virtual Private Network\) and shared drives \(German only\)](#)
- ↳ [Granting temporary access to TUD WLAN](#)
- ↳ [Neuroimaging Center Python Pipeline \(NICePipe\)](#)
- ↳ [NIC XNAT Database](#)
- ↳ [PaMS](#)
- ↳ [Redcap Study Database at the ZIH](#)

### Open Science Resources


- ↳ [Open Access Repository and Archive \(OpARA\)](#)
- > [Open Science Initiative der Fakultät Psychologie der TUD \(OSIP\)](#)
- ↳ [Open Science Framework \(OSF\)](#)

<https://nic-tud.de/research/ressourcen/data-management-manual/>

## NIC Data-Management Manual

### Procedure for study staff:

Click on the individual steps for detailed instructions.




Please send feedback to [mrphysik@msx.tu-dresden.de](mailto:mrphysik@msx.tu-dresden.de)

# Outline

---





- MRI data flow
- Quality flags for sequences
- Upload of additional data
- DICOM→NIFTI Converter
- Download data
- User groups (Owner, Member, Collaborator)
- Outlook: Automatic image processing software: NiCePype

# MRI Data Flow



# MRI Data Flow: PaMS registration

Neuroimaging Zentrum — Probandendatenbank

Christian Löschen    

Start Probandeninfos **Untersuchungen** Befundungen Implantate Anrufprotokoll Unbefundet Befundet Rekrutieren Studien Berechtigungen ORSEE Scanner-Reboot

Spur: Implantate (Test Basicassessment) • Anrufprotokoll (Test Basicassessment) • Unbefundet • Befundet • Rekrutieren • Studien • Berechtigungen • Start • Probandeninfos (Test Basicassessment) • Studien (Test Basicassessment)

**Proband #000844**  
Test Basicassessment (17)

**Neue Untersuchung**

Untersuchungsart  
Magnetresonanztomographie (MRT)

Scanner: 1055 Untersuchungsraum: A/11

Gewicht [kg]:  Körpergröße [m]:

Studiennummer:  Probandencod:





Studienkommentar:

Erhebungsnummer:  Untersuchungsdatum: 06.02.2020

1. Untersucher: Löschen, Christian 2. Untersucher: Aguilera, Alvaro

Untersuchungsstatus:  Teilnahmestatus:

**Untersuchungen**

Untersuchungs-ID	Studie	Participant ID	Erhebungsnummer	Datum	Untersucher #1	Status
1	10101	99999	1	06.02.2020 00:00	Christian Löschen	   

Export Registration to scanner

3	1	2	
1	Surname	Study number	10101
2	Given name	Participant code	99999
3	Subject no.	PaMS-ID (6 digits)	000844

# Automatic import to the scanner console

**Patientenregistrierung**

**PATIENT**

**1** Nachname 10101

**2** Vorname 99999

Titel 1

**3** Patientennr. 000844

Geburtsdatum 25.06.1971 [dd.MM.yyyy]

Geschlecht  männl.  weiblich  Andere

Alter 48 Jahre

Größe 1,83 [m]

Gewicht 72,00 [kg]

Bemerkungen Project:10101;Subject:99999;Sex

**METHODE**

Auftragsnummer 14915

Anforderungsnr.

Erforderliche Methode(n)

Studie Kopf TRR265

Studienkommentar

Patientenlage Kopf voran - Rückenlage

Institut Systemische Neurowissenschaft

**HOSPITAL**

Überweisender	1	Surname	Study number	10101
Anfordernder	2	Given name	Participant code	99999
Überweisung	3	Subject no.	PaMS-ID (6 digits)	000844

Vormerken Untersuchen Suchen Abbrechen Hilfe

ISO\_IR 100

# XNAT Levels

TECHNISCHE UNIVERSITÄT DRESDEN

## Project Overview

XNAT\_ZIH currently contains 186 projects, 3176 subjects, and 3635 imaging sessions.

Projects Subjects MR PET CT

ID:  Name:

Keywords:  Investigator: (SELECT)

**Projects**

- Aging and dopamine modulation of complementary control processes  
Project ID: 00034  
You are an owner for this project.
- NICDB Testprojekt  
Project ID: 01333 Pt: Dirk Müller  
You are a member for this project.
- Neuromodulation of structural and functional changes by homeostatic signals after smoking cessation  
Project ID: 00032  
You are an owner for this project.
- Neuronale Plastizität im Rahmen einer kognitiv-verhaltens-therapeutischen Intervention bei Patienten mit Panikstörung  
Project ID: 00201  
You are an owner for this project.
- Serotonergic modulation of meta-control parameters  
Project ID: 00030  
You are an owner for this project.

## NICDB Testprojekt Participant Overview

**Details**

ID: 01333  
PI: Müller, Dirk

**Actions**

- Add
- Upload Images

**Subjects** ✕

<< first < prev 1 next > last >> 200 1 of 1 Pgs (16 Rows)

Subject	M/F	Hand
006217	U	
20190215	U	
20190624	U	
20190625	U	
20191203	U	
20191213	U	

PROJECT: 01333 > SUBJECT: 20190624 > 006217

## MR Session: 006217 Session Overview

**Details** Projects

Accession #: XNAT\_ZIH2\_E01301 Subject: 20190624  
Date Added: 11/12/2020 13:28:31 (rbaensch) Gender:  
Date: 06/24/2019 Handedness:  
Time: 11:22:00 Age: --  
Scanner Name: MRC35213  
Scanner Type: SIEMENS TrioTim  
Acquisition Site: NIC

**Actions**

- Edit
- Build
- View
- Download
- Email
- Manage Files
- View Images
- Upload Additional Files

**Scans**

Bulk Actions: Download

Scan	Type	Series Desc	Usability	Files	Note
<input type="checkbox"/> 1	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	usable	2.8 MB in 3 files	
<input type="checkbox"/> 2	boldA_restfb_01_32traa-6117-20202000-0987-33	boldA_restfb_01_32traa-6117-20202000-0987-33	usable	30.7 MB in 20 files	
<input type="checkbox"/> 3	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	usable	1.5 MB in 1 files	
<input type="checkbox"/> 4	boldA_restfb_01_32traa-6117-20202000-0987-33	boldA_restfb_01_32traa-6117-20202000-0987-33	usable	10.7 MB in 7 files	

Total: 45.8 MB in 31 files

# Pre-Archive

- Data are automatically transferred from the MRI to the Pre-Archive
- Images are *exclusively* evaluated by the MTRAs for completeness
- Complete data sets are sent to the project archive and marked as *unchecked*
- Study personal ought to work only with the transferred data in the respective project archive!


The screenshot displays a web application interface for managing MRI sessions. The main area is a table titled 'Sessions' with columns for Project, Subject, Session, Scan Date, Upload Date, and Status. A single row is visible with the following data: Project: 00034, Subject: 202060, Session: 16358\_1, Scan Date: Tue Dec 01 2020..., Upload Date: Tue Dec 01 2020 12..., Status: Ready. Above the table are buttons for 'Refresh', 'Check', 'Uncheck', 'Check All', and 'Uncheck All'. To the right of the table is a 'Details of Highlighted Session(s)' panel, which currently shows 'No sessions highlighted.' Below this panel are 'Actions' (Archive, Review and Archive, Change Projects, Delete, Rebuild) and 'Filters' (Projects, Upload Date, Scan Date, Status, Subject/Session). A large red watermark 'MTRAs ONLY!' is overlaid on the table area.

Project	Subject	Session	Scan Date	Upload Date	Status
00034	202060	16358_1	Tue Dec 01 2020...	Tue Dec 01 2020 12...	Ready



# Step 1: Quality flags - Handling incorrect or unwanted data

- Scans for positioning the subjects
  - ➔ AAScout
  - ➔ Localizer
- Data from unwanted reconstructions
  - ➔ T1 thick film
  - ➔ DTI
- Scan aborts at the request of the subjects and repeated measurements



# Step 1: Archive: Set Quality Flag

Details Projects

**Accession #:** XNAT\_ZIH2\_E01590      **Subject:** 20201203-35213  
**Date Added:** 12/04/2020 11:27:56 (mame720c)      **Gender:**  
**Date:** 12/03/2020      **Handedness:**  
**Time:** 07:06:54      **Age:** --  
**Operator:** M?ller Dirk  
**Scanner Name:** MRC35213  
**Scanner Type:** SIEMENS TrioTim  
**Acquisition Site:** NIC

**Actions**

- Edit
- Build
- View ▶
- Download ▶
- Email
- Manage Files
- View Images
- Delete

## Scans

Bulk Actions: Download

<input type="checkbox"/>	Scan	Type	Series Desc	Usability	Files	Note
<input type="checkbox"/>	1	ep2d_bold_rsMRT_Phantom	ep2d_bold_rsMRT_Phantom	UNCHECKED	<a href="#">4.1 MB in 11 files</a>	
<input type="checkbox"/>	2	ep2d_bold_rsMRT_Phantom	ep2d_bold_rsMRT_Phantom	UNCHECKED	<a href="#">55.8 MB in 150 files</a>	
<input type="checkbox"/>	3	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	UNCHECKED	<a href="#">1.5 MB in 1 files</a>	
<input type="checkbox"/>	4	boldA_restfb_01_32traa-6117-20202000-0987-33	boldA_restfb_01_32traa-6117-20202000-0987-33	UNCHECKED	<a href="#">384.0 MB in 250 files</a>	
<input type="checkbox"/>	5	boldA_TRR-rest_1-1_32t-a-1-6-1-8-24242400-38-869_SBRef	boldA_TRR-rest_1-1_32t-a-1-6-1-8-24242400-38-869_SBRef	UNCHECKED	<a href="#">1.1 MB in 1 files</a>	
<input type="checkbox"/>	6	boldA_TRR-rest_1-1_32t-a-1-6-1-8-24242400-38-869	boldA_TRR-rest_1-1_32t-a-1-6-1-8-24242400-38-869	UNCHECKED	<a href="#">263.1 MB in 250 files</a>	

Total: [709.5 MB in 663 files](#)

# Step 1: Archive: Set Quality Flag

## Scans

Add Scan

Scan	Type ?	Quality	Note	
1	AAHScout_32	usable		15.2 MB in 128 files
2	AAHScout_32_MPR	(SELECT) UNCHECKED		491.9 KB in 4 files
3	localizerA_standard_1_32-m	usable		1.1 MB in 7 files
4	T1wB_standard_1_32-s-a-1	questionable unusable		65.0 MB in 240 files
5	T1wD_vibeFlip15_1_32-t-r-2	usable		15.2 MB in 60 files
6	T1wD_vibeFlip15_2_32-t-r-2	usable		15.2 MB in 60 files
7	T1wD_vibeFlip15_2_32-t-r-2	usable		15.2 MB in 60 files

Total: 127.4 MB in 559 files


## Step 1: Archive: Set Quality Flag

---

- **Usable:** data remain in archive and can be used for further processing
- **Questionable:** data remain in archive but should be reviewed before further processing
- **Unusable:** data will be removed from the project archive but can be recovered if necessary

## Step 2: Upload of additional data

- Paradigm logging files
- Pulse oximetry
- Electrodermal activity
- ...



# Step 2: Upload of additional data

[PROJECT: 01333](#) > [SUBJECT: 20200527](#) > [15133\\_1](#)

MR Session: 15133\_1

Details		Projects	
<b>Accession #:</b>	XNAT_ZIH02_E02675	<b>Subject:</b>	<a href="#">20200527</a>
<b>Date Added:</b>	05/27/2020 11:50:10 (admin)	<b>Gender:</b>	
<b>Date:</b>	05/27/2020	<b>Handedness:</b>	
<b>Time:</b>	11:39:10	<b>Age:</b>	--
<b>Scanner:</b>	MRC35213 SIEMENS TrioTim		
<b>Acquisition Site:</b>	NIC		

**Actions**

- Edit
- View ▶
- Download ▶
- Email
- Manage Files**
- View Images
- Delete

## Scans

Scan	Type	Series Desc	Usability	Files	Note
⊕ 2	AAScout	AAScout	usable	23.0 MB in 256 files	
⊕ 4	ep2d_bold_rsfMRT_Phantom	ep2d_bold_rsfMRT_Phantom	usable	2.3 MB in 6 files	
				<b>Total:</b> 25.2 MB in 262 files	

# Step 2: Upload of additional data

The screenshot shows a 'File Manager' window with a tree view of folders. The 'scans' folder is expanded, showing sub-folders: '2 AAScout', 'DICOM 256 files, 22.96 MB DICOM RAW', 'SNAPSHOTS 2 files, 161 KB GIF SNAPSHOTS', and '4 ep2d\_bold\_rsMRT\_Phantom'. The '4 ep2d\_bold\_rsMRT\_Phantom' folder is expanded, showing 'DICOM 6 files, 2.25 MB DICOM RAW' and 'SNAPSHOTS 2 files, 99 KB GIF SNAPSHOTS'. An 'Add Folder' dialog box is open on the right, with 'Level' set to 'scans', 'Item' set to '2', and 'Folder' set to 'Logfiles'. The 'Format Information' section has empty fields for 'Format', 'Content', and 'Tags'. At the bottom right of the dialog are 'Create' and 'Close' buttons. At the bottom of the File Manager window, there are buttons for 'Add Folder', 'Upload Files', 'Update File Data', a 'zip' dropdown, 'Download', and 'Close'. The 'Add Folder' button is circled in red.

Upload additional data such as paradigm logging files or physiological data

# Step 2: Upload of additional data

The image shows a web interface with two main panels: 'Upload File' and 'File Manager'.

**Upload File Panel:**

- Destination:** Level: scans, Item: 2, Folder: Logfiles.
- File Information:** Fields for Rename, Format, Content, and Tags. An Overwrite checkbox is present.
- Search:** A search bar containing 'Logfile\_sj\_9999.txt'.
- Buttons:** 'Upload' and 'Close'.

**File Manager Panel:**

- Tree view showing folders: Resources, paradigm (1 files, 34 bytes), scans, 2 AAScout, Logfiles (1 files, 0 bytes), DICOM (256 files, 22.96 MB DICOM RAW), SNAPSHOTS (2 files, 161 KB GIF SNAPSHOTS), 4 ep2d\_bold\_rsfMRT\_Phantom, log\_files (1 files, 0 bytes), SNAPSHOTS (2 files, 99 KB GIF SNAPSHOTS), paradigm (1 files, 34 bytes gif cross), and DICOM (6 files, 2.25 MB DICOM RAW).
- The file 'Logfile\_sj\_9999.txt' (0 bytes) is highlighted with a red circle.
- Bottom navigation buttons: 'Add Folder', 'Upload Files' (circled in red), 'Update File Data', 'zip', 'Download', and 'Close'.



# Step 3: XNAT pipelines


Details Access Manage **Pipelines**

ID: 00034

Edit Details Delete Manage Custom Variables

**Actions**


- Add
- Upload Images
- View Prearchive
- Scan Type Cleanup
- Add to Favorites
- Download XML
- Download Images
- Manage Files



**Project owners only!**


# Step 3: XNAT pipelines

Aging and dopamine modulation of complementary control processes



The screenshot shows a navigation bar with four tabs: 'Details', 'Access', 'Manage', and 'Pipelines'. The 'Pipelines' tab is highlighted in blue and circled in red. Below the navigation bar, the ID '00034' is displayed. Three buttons are visible: 'Edit Details', 'Delete', and 'Manage Custom Variables'.

Aging and dopamine modulation of complementary control processes



The screenshot shows the 'Pipelines' tab selected and highlighted in blue. Below the navigation bar, the text 'Pipelines for 00034' is displayed. A table with the following columns is shown: 'Applies To', 'Generates', 'Name', and 'Description'. The table contains the text 'No records found.' Below the table, the button 'Add More Pipelines' is circled in red.

**Project owners only!**

# Step 3: XNAT pipelines

Aging and dopamine modulation of complementary control processes

Details Access Manage Pipelines

*Additional Pipelines for 00034*

		Applies To	Generates	Name	Description
<a href="#">Add</a>	<a href="#">Details</a>	All Datatypes		dicom2nifti.xml	New DICOM to NIFTI image data converter (dcm2niix)

Show Project Pipelines

**Project owners only!**

# Step 3: Convert DICOMs to NIFTI


[PROJECT: 01303](#) > [SUBJECT: 20200221](#) > [15012\\_1](#)

MR Session: 15012\_1

Details		Projects	
<b>Accession #:</b>	XNAT_E00755	<b>Subject:</b>	20200221
<b>Date Added:</b>	02/21/2020 06:42:49 (admin)	<b>Gender:</b>	
<b>Date:</b>	02/21/2020	<b>Handedness:</b>	
<b>Time:</b>	06:29:18	<b>Age:</b>	--
<b>Operator:</b>	Geißler		
<b>Scanner Name:</b>	MRC35213		
<b>Scanner Type:</b>	SIEMENS TrioTim		
<b>Acquisition Site:</b>	NIC		

Actions
Edit
View ▶
Download ▶
Email
<b>Build</b>
Manage Files
View Images
Delete

# Step 3: Convert DICOMs to NIFTI




# Step 3: Convert DICOMs to NIFTI

New DICOM to NIFTI image data converter (dcm2niix) for 15012\_1

Name	Value
scanIDs	<input checked="" type="checkbox"/> 1 <input checked="" type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input checked="" type="checkbox"/> 5
sessionID	<input type="text" value="XNAT_E00755"/>
sessionName	<input type="text" value="15012_1"/>
projectID	<input type="text" value="01303"/>
subjectID	<input type="text" value="XNAT_S00649"/>
compress	<input type="text" value="n"/>
notify	<input type="text" value="0"/>

- Run pipeline in-process?**
- Suppress workflow entries for pipeline execution?**

# Step 3: Convert DICOMs to NIFTI



# Download Data

MR Session: 006217

**Details** | Projects

**Accession #:** XNAT\_ZIH2\_E01301      **Subject:** 20190624  
**Date Added:** 11/12/2020 13:28:31 (rbaensch)      **Gender:**  
**Date:** 06/24/2019      **Handedness:**  
**Time:** 11:22:00      **Age:** --  
**Scanner Name:** MRC35213  
**Scanner Type:** SIEMENS TrioTim  
**Acquisition Site:** NIC

**Actions**


- Edit
- Build
- View ▶
- Download ▶
- Email
- Manage Files**
- View Images
- Upload Additional Files

## Scans

**Bulk Actions:** Download


<input type="checkbox"/>	Scan	Type	Series Desc	Usability	Files
<input type="checkbox"/>	1	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	usable	<a href="#">2.8 MB in 3 files</a>
<input type="checkbox"/>	2	boldA_restfb_01_32traa-6117-20202000-0987-33	boldA_restfb_01_32traa-6117-20202000-0987-33	usable	<a href="#">30.7 MB in 20 files</a>
<input type="checkbox"/>	3	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	boldA_restfb_01_32traa-6117-20202000-0987-33_SBRef	usable	<a href="#">1.5 MB in 1 files</a>
<input type="checkbox"/>	4	boldA_restfb_01_32traa-6117-20202000-0987-33	boldA_restfb_01_32traa-6117-20202000-0987-33	usable	<a href="#">10.7 MB in 7 files</a>

Total: 45.8 MB in 31 files





# Download Data



# User rights management

## Understanding User Roles and Permissions

By default, XNAT defines three common project roles: **Owners**, **Members**, and **Collaborators**. As a project owner, you have the ability to add specific users to your project within these roles. These users will then have the permissions associated with their defined role (which overrules the default accessibility permissions).

Role/Activity	Project Owners	Project Members	Project Collaborators
Create Data	<b>C</b>	<b>C</b>	
Read/Download Data	<b>R</b>	<b>R</b>	<b>R</b>
Update Data	<b>U</b>	<b>U</b>	
Delete Data	<b>D</b>		

*The "CRUD" model of permissions*

- **Project Owners:** If you define additional project owners, they will have all of the permissions on your project that you do. They can read, insert, modify, and delete anything (and everything) associated with your project. They can also add additional users to your project and modify the data types associated with your project.
- **Project Members:** Members have the ability to manage the data in your project. They can read, insert, and modify subjects and experiments in your project. They cannot modify the project users and data types.
- **Project Collaborators:** Collaborators have read-only access on all of the data in your project. They cannot insert or modify data owned by your project. They can download your data and use it within their projects.

# Outlook: Automatic MRI Data Preprocessing: NiCePype

Currently under development for the TRR265

Specify required data and evaluation pipelines once only

Last login: 12/09/2020 17:49:48

Browse ▾ New ▾ Upload ▾ Administer ▾ Tools ▾ Help ▾

**TECHNISCHE UNIVERSITÄT DRESDEN**

XNAT\_ZIH currently contains 186 projects, 3177 subjects, and 3636 imaging sessions.

Projects Subjects MR PET CT

Automatic download of imaging data

```
[10101|sst_1_2_61_3|MRC35396]
dir_name=mannheim
subject_regex=RecoTest
fmri_magnitude_regex=bold\S*SST
fmri_magnitude_frames=200
fieldmap_magnitude_regex=fmapA\S*
fieldmap_phase_regex=fmapA\S*
```

## NiCePype

Neuroimaging Center Python Pipelines

Overview Data Processing-Jobs Logout

Project Overview: 1301 (1301)

Found a total of 0 data processing jobs for this project in the database, whereof 0

Project	1301 (1301)
Department	imphys


Processed data will be automatically uploaded, including QA

# Outlook: Automatic MRI Data Preprocessing: NiCePype

Currently under

Last login: 12/09/2020 17:49:48


Browse ▾ New ▾ Upload



TECHNISCHE UNIVERSITÄT DRESDEN

XNAT\_ZIH currently contains 186 projects

Projects Subjects MR



File Manager

- Resources
  - sst\_1\_2\_60\_3:dresden 198 files, 4.04 GB NIFTI
    - nifti\_dcm2nii NIFTI
    - preprocessing TXT
    - qa\_report PNG
  - scans
    - 1 AAHScout\_32
      - DICOM 128 files, 15.22 MB DICOM RAW
      - SNAPSHOTS 2 files, 315 KB GIF SNAPSHOTS
    - 2 AAHScout\_32\_MPR
      - DICOM 4 files, 492 KB DICOM RAW
      - SNAPSHOTS 2 files, 25 KB GIF SNAPSHOTS
    - 3 localizerB\_standard\_1\_32-m-a-1-1-1-8-191
      - DICOM 7 files, 1.10 MB DICOM RAW
      - SNAPSHOTS 2 files, 26 KB GIF SNAPSHOTS
    - 4 T1wA\_biobank\_1-1\_32-s-a-1-1-2-88-10101
      - DICOM 208 files, 41.62 MB DICOM RAW
      - SNAPSHOTS 2 files, 748 KB GIF SNAPSHOTS
    - 5 boldA\_TRR-rest\_1-1\_32-t-a-1-6-1-8-24242400-38-869\_SBRef
      - DICOM 1 files, 1.07 MB DICOM RAW
      - SNAPSHOTS 2 files, 78 KB GIF SNAPSHOTS

File Manager

- Resources
  - qa\_report\_sst\_1\_2\_61\_3\_mannheim 1 files, 1.34 GB
    - qa\_report\_sst\_1\_2\_61\_3\_mannheim.zip 1.34 GB
  - qa\_report\_sst\_1\_2\_61\_3\_berlin 1 files, 706.75 MB

Specify required data and evaluation pipelines once only

```
[sst_1_2_61_3|MRC35396]  
name=mannheim  
t_regex=RecoTest  
magnitude_regex=bold\S*SST  
magnitude_frames=200  
map_magnitude_regex=fmapA\S*  
map_phase_regex=fmapA\S*
```

Processed  
automatic  
including (

Project in the database, whereof	
	1301 (1301)
	imphys

# Thank you

## Ressources

Contact: [sfb940-inf@groups.tu-dresden.de](mailto:sfb940-inf@groups.tu-dresden.de)  
[mrphysik@msx.tu-dresden.de](mailto:mrphysik@msx.tu-dresden.de)

NIC homepage:  
<https://nic-tud.de/research/ressourcen/>

XNAT webpage:  
<https://xnat.zih.tu-dresden.de>

Select „TU Dresden“ as Login and use your  
ZIH Login

