



## Terms and conditions for the NMR sample submission

- Generally, only cleaned and intact NMR sample tubes are allowed. The sample will not be measured when the NMR sample tube is whether dirty and/ or damaged.
- Samples for which the required effort exceeds the standard conditions or measuring times (2D NMR, hetero nuclei) have to be announced in advance. Please ask the corresponding NMR facility manager.
- Only designated NMR tubes are allowed (5 mm diameter). The minimum lengths must be more than 15 cm.
- NMR caps which are porous, cracked or otherwise damaged (dirty) are not allowed.
- NMR sample tubes must be wiped off (with e.g. ethanol) – no grease on the outer glass-wall.
- Parafilm or teflon stripes must be adjusted properly. Please no remains or rims.
- Samples with broken glassware (broken pipette) inside the tube are not allowed.
- NMR samples which are potentially under pressure are not allowed. In case, please use a J Young tube. This must be announced in advance.
- The NMR sample sheet must be filled out properly (affiliation, deuterated solvent, experiment). Otherwise, the sample will not be measured.
- Please pick-up your sample in time after you have received the NMR data.
- By handing out the NMR sample sheet, you accept the aforementioned criteria and the accounting.